

# PROTOCOL: Final Sample Resuspension (FSR)

## Purpose

Prepare samples for mass spec analysis following the final desalt on the AssayMAP BRAVO. The addition of heavy synthetic peptides allow for QC during analysis, quantify phospho-sites in endogenous samples, and compare phospho-mark abundance across sample sets.

This protocol is to be followed the day samples are ready to be analyzed on the Q Exactive HF

## Preparation

1. Thaw Heavy Synthetic Peptide Mixture {FSR-M04} on ice
2. Remove dried, desalted, eluate from the freezer location. Bring to room temperature.

## Materials

- Acetonitrile (ACN), EMD Millipore, Cat. No. AX0156-1 {FSR-M01}
- Formic Acid (FA), Sigma-Aldrich, Cat. No. S6302-50ML-F {FSR-M02}
- HPLC-grade water, JT Baker, Cat. No. 4218-03 {FSR-M03}
- Concentrated Heavy Synthetic Peptide Mixture {FSR-M04}
- Micronic Split-Septa TPE Push Caps, Micronic, Cat. No. 1775-3026 {FSR-M05}
- 1L Glass Graduated Cylinder {FSR-M06}
- 1L Glass Reagents Bottle {FSR-M07}
- 40 mL Amber Vial, Solid Screw Top, PTFE Liner {FSR-M08} [SUPLECO, 27182]
- 1.7 mL Microcentrifuge Tubes {FSR-M09} [Axygen, MCT-175-L-C]
- 0.5 mL Microcentrifuge Tubes {FSR-M10} [Axygen, MCT-060-L-C]

## Assets

## Reagent Mixes

ID	Name	Step	Composition	Stock Volume	Use
FSRMIX01	3%ACN/5%FA	FSR	3% {FSR-M01} and 5% FA {FSR-M02} in {FSR-M03}	40 mL	Balance reagent in creation of heavy synthetic peptide mixture
FSR MIX02	1:50 Dilution of Concentrated Heavy Synthetic Phosphopeptide Mix {FSR-M04}	FSR	1:50 dilution of concentrated heavy phosphopeptides in 3%ACN/5%FA {FSRBUF01}	22uL at 50x concentration	Mixture of heavy synthetic peptides used as an internal standard in the analysis of samples

## Preparation of Mixture Reagents:

FSRMIX01- 3% ACN/ 5% FA :

- Add 36.8 mL of HPLC-grade water {FSR-M03} to 40 mL amber vial with screw top{FSR-M10}
- Add 2 mL of formic acid {FSR-M02} to vial
- Add 1.2 mL of acetonitrile {FSR-M01} to vial
- Vortex to mix.

FSRMIX02- 1:50 Dilution of Heavy-tagged Phosphopeptide Standard Mix

- Add 1078uL of 3%ACN/5%FA {FSRMIX01} to 22uL of Concentrated Heavy-tagged Phosphopeptide Standard Mix {FSR-M04}.

## Procedure

1. Deliver 10uL of FSRMIX02 to each sample.
2. Cover the individual sample vials with split-septa TPE push caps {FSR-M05}.
3. Spin down the plate.
4. Keep plate at 4°C until it is ready to be placed on the the deck of the HPLC system.