

Study 9-1 SOP (Phase II) for NCI CPTAC Consortium-Wide Multiple Reaction Monitoring (MRM) Experiment 9-Point Calibration Curve of 123 target peptides

Experimental Design and Statistics Verification Studies Working Group

Overview

Study 9-1, the 9 point calibration curve will entail monitoring 123 signature peptides and corresponding stable-isotope labeled internal standard (SIS) peptides (~750 total transitions) that will be quantitatively assayed by time-scheduled LC-MRM-MS against a background of MARS-14 depleted human K₂EDTA plasma (0.5 µg/µL).

Aims: The study described in this SOP is designed to accomplish the following:

1. To generate a 9 point standard curve for 123 potentially cancer relevant peptide targets spanning 1 amol/µL to 100 fmol/µL in a depleted digested plasma background (500 ng/µL) that will be analyzed by LC-MRM-MS on four different triple quadrupole platforms at volunteer CPTAC sites.
2. To enable calculations of LOD, LOQ, accuracy and precision for all peptides from the calibration curve.
3. A set of 6 blinded samples is provided that mimic real biomarker samples.
4. To evaluate the level of endogenous and interfering signals for all peptides by replicate measurements of blank plasma samples.

The number of transitions that will be targeted in this method (~750) is a critical challenge and will require the use of retention time scheduled LC-MRM-MS. To achieve success, the following protocol is designed to assess the system suitability and retention time stability of the LC-MRM-MS instrument, as well as design a single MRM-MS method empirically. The protocol outlined herein should be followed as strictly as possible and all deviations from this protocol must be outlined in detail on the last sheet. The experimental design is outlined below.

Experimental Outline

1. Run System Suitability Sample 5 times to condition column and assess LC-MRM-MS performance
2. Using Skyline, create 6 MRM-MS transition lists for the heavy target peptides only (maximum of 60 transitions in each unscheduled method)
3. Analyze the heavy-only peptide sample using the 6 different MRM-MS methods
 - i. Import the 6 raw files into Skyline
 - ii. Obtain retention times for each heavy target peptide
 - iii. Export a scheduled MRM-MS transition list (or MRM-MS method) with a 4 min RT window
 - iv. Run the scheduled LC-MRM-MS method twice using the sample containing heavy peptides spiked into depleted, digested plasma.
 - v. Import the 2 raw files into Skyline
 - vi. Add in the light (unlabeled) peptides into the Skyline document
 - vii. Export a new scheduled LC-MRM-MS method with all ~750 transitions and a 2 min RT window
4. Execute calibration curve and blinded samples according to specified run order.
5. Import data into Skyline and integrate.
 - i. Import System Suitability Sample files into System Suitability Skyline file and integrate. Ensure peak area CVs and retention times are within specified error limits.
 - ii. Import Calibration Curve files into Study 9-1 Skyline file and integrate.
6. Export reports from each Skyline file and upload to the NIST ftp server (see separate document).

Materials and Reagents

A. System Suitability Sample (Study 9-1-SSS)

- a. Digest of 6 equimolar proteins (Michrom Bioresources, #PTD/00001/63)
 - i. Four 10 μL aliquots, 1 pmol/ μL , supplied in 30% acetonitrile/0.1 % formic acid in water

B. Heavy-Only IS Peptide Mixture (in 25 fmol/ μL 6 protein mix)

- a. 123 heavy IS peptides at 800 fmol/ μL each
 - i. Two 10 μL aliquots, supplied in 30% acetonitrile/0.1% formic acid

C. QC Sample (Sample 9-1-QC)

- a. Equimolar mixture of the 123 unlabeled and 123 labeled synthetic peptides in 25 fmol/ μL 6 Protein mix matrix
 - i. Four 25 μL aliquot supplied at 10 fmol/ μL each peptide in 0.1% formic acid in water. Background of 25 fmol/ μL of 6 protein mix added to minimize adsorption/loss of hydrophobic peptides.

D. Digested human plasma spiked with labeled IS peptides only, (Sample 9-1-A)

- a. Depleted plasma diluted to a final concentration after digestion of approximately 0.5 $\mu\text{g}/\mu\text{L}$
 - i. 123 labeled IS peptides spiked at a concentration of 10 fmol/ μL (Sample 9-1-A)
 - ii. Ten 25 μL aliquots supplied in 0.1 % formic acid in water (after desalting by SPE, note: desalting was already performed at Vanderbilt University)

E. Digested human depleted plasma spiked with 123 unlabeled synthetic peptides and 123 labeled IS peptides (Samples 9-1-B to 9-1-J)

- a. Depleted plasma diluted to a final concentration after digestion of ~ 0.5 $\mu\text{g}/\mu\text{L}$
 - i. 123 labeled IS peptides spiked at a concentration of 10 fmol/ μL
 - ii. Four 25 μL aliquot of each spike level supplied for 4 singlicate curves
 - iii. 123 unlabeled synthetic peptides are spiked in at the following approximate concentrations:

| Sample (Study 9-1) | Spiked [$^{12}\text{C}/^{14}\text{N}$] peptide Concentration (fmol/ μL) |
|-----------------------|---|
| 9-1-J | 100.00 |
| 9-1-I | 23.71 |
| 9-1-H | 5.62 |
| 9-1-G | 1.33 |
| 9-1-F | 0.316 |
| 9-1-E | 0.075 |
| 9-1-D | 0.018 |
| 9-1-C | 0.004 |
| 9-1-B | 0.001 |

Important note: each of the 4 different singlicate calibration curves will be completed before the next replicate curve will be started: Specific naming including the first blank (A sample) is required, for details also see run order below:

9-1-A1 to 9-1-J1; 9-1-A2 to 9-1-J2; 9-1-A3 to 9-1-J3; 9-1-A4 to 9-1-J4.

F. Six Blinded Samples to be run in four replicates in between each singlicate calibration curve (Total: Samples 9-1-01Blinded to 9-1-24Blinded)

Four replicates of 6 Blinded samples (each provided in 25 aliquot):

first replicate: 9-1-01Blinded to 9-1-06Blinded
second replicate: 9-1-07Blinded to 9-1-12Blinded
third replicate: 9-1-13Blinded to 9-1-18Blinded
fourth replicate: 9-1-19Blinded to 9-1-24Blinded

mixed run order for the Blinded Samples within each replicate set

G. "Wash Samples", to avoid carryover between blinded samples, the run order defines running a "wash" gradient in between blinded samples. "Wash Samples" will be 10 fmol/uL of the 6 protein mix (Michrom sample). Samples will be diluted from provided 6 protein mix / Michrom sample at stock concentration of 1 pmol/ μ L).

H. HPLC solvents

- a. Acetonitrile, HPLC grade
- b. Water, HPLC grade
- c. Formic Acid

I. PicoFrit columns

- a. 75 μ m ID with 10 μ m ID tip
- b. For ABSciex and Thermo operators: **Prepacked New Objective Reprosil Columns** that were sent to all laboratories from Vanderbilt University (in case the operator runs out of prepacked columns please contact Susan Abbatiello or Birgit Schilling)
- c. For Waters UPLC operators: a Waters specific column will be used (1.7 μ m BEH130 C18, 75 μ m X 150mm column)

J. Loop for Eksigent LC: 1 μ L PEEKsil loop provided in sample kit, 100 μ m ID PEEKsil; newer design 1 μ L loops are provided for Eksigent Ultra systems.

K. Fused silica

- a. 20 μ m ID between gradient pump and autosampler (for direct inject configuration)
- b. 20 μ m ID between autosampler and PicoFrit column (keep tubing as short as possible)

L. High voltage contact

- a. For liquid-liquid junction, Upchurch part P-888.

M. Autosampler vials

- a. Polypropylene autosampler vials with conical insert, 250 μ L maximum volume

N. Autosampler loop

- a. 1 μ L loop for direct inject set-up (fused silica or PEEKsil)

O. Column heater

- a. If the system is equipped with a column heater, please set to 35 degrees Celsius.

Required Software

- A. **Skyline Daily Software**, most recent version (please note version used)
- B. **Platform-specific Skyline document** (*.sky), provided by CPTAC VWG
- C. **Microsoft Excel** (or similar)
- D. **AuDIT Software** (available at <http://www.broadinstitute.org/cancer/software/genepattern/modules/AuDIT.html>)

Overview

In Main Study 9-1 (9-point calibration curve), synthetic [$^{12}\text{C}/^{14}\text{N}$] and [$^{13}\text{C}/^{15}\text{N}$] signature peptides will be spiked into MARS-14 depleted, digested human K_2EDTA plasma and analyzed by scheduled LC-MRM-MS at all sites. All sample preparation will be performed at Vanderbilt University prior to distribution of the sample kits. Results from Study 9-1 from all CPTAC sites, will provide LOD/LOQ values for each peptide (this design contains minimal sample handling, and no target protein digestion and eliminates key factors of experimental variability). The Study 9-1 scheduled LC-MRM-MS samples will consist of a 9-point standard curve ranging in concentration from 100 fmol/ μL to 1 amol/ μL , with all concentrations in the presence of 0.5 $\mu\text{g}/\mu\text{L}$ depleted digested plasma. Synthetic signature internal standard (IS) peptides uniformly labeled with a $^{13}\text{C}/^{15}\text{N}$ amino acid will be spiked into all plasma samples at a constant concentration of 10 fmol/ μL .

The following paragraphs provide a description of the experimental procedure, with each major step separated into its own paragraph and heading. Paragraphs are followed by outlined “checklists” that have matching headings for reference.

1. System Suitability Sample, Column Conditioning and Instrument Performance

System Suitability Sample (SSS) runs using the Michrom 6 protein mix (based on Study 9S) will be interspersed into Study 9-1 (to guarantee system suitability and performance and, in particular, monitor peak area stability and potential RT drift). The SSS will first be analyzed in an unscheduled LC-MRM-MS method (Study 9S) to condition the column and assess the performance of the LC-MS instrument platform. Upon obtaining the specified results for retention time and peak area, sites may proceed to the generation of the scheduled LC-MRM-MS runs for Study 9.

The SSS will also be run every 6-8 runs in the sample queue, in order to track instrument performance. Data will be analyzed in Skyline and will be submitted in a separate “Michrom Study 9S” file.

2. Generation of Scheduled LC-MRM-MS Method for Study 9 Peptides

In preparation for the scheduled calibration curve LC-MRM-MS runs, in which 738 transitions will be monitored in one run, participating sites will monitor all heavy synthetic labeled peptides with about 369 transitions (123 peptides with 3 transitions each) in six unscheduled LC-MRM-MS runs (~60 transitions per run), with each run represented by a separate LC-MRM-MS method. Retention times for scheduling will be determined empirically using these six runs and will be verified by analysis with a single scheduled LC-MRM-MS run prior to analyzing the samples from which the calibration curve will be generated. All method building and data analysis will be performed using Skyline.

3. Data Analysis

The experimentally determined molar concentration of the spiked peptide or protein will be calculated and compared to its theoretical value for accuracy. Operators will import all acquired data files into Skyline and check and if necessary adjust peak integration. Check "Integrate All" (Menu, Settings, Integrate all), so that heavy and light transitions will be "integrated together".

While checking peak integration open RT replicate views (Menu, View, Retention Time, Replicate Comparison), and Peak Area replicate views (Menu, View, Peak Areas, Replicate Comparison) to visually help confirm proper peak integration. Specifically, check for interferences and RT drift problems during scheduling. Use Skyline Custom Annotation features to annotate any observations or notes for Study Statisticians to consider (Menu, View, Results Grid; and annotate observations on the precursor or transition level). Finally use "AuDIT" to analyze data for interferences before submitting data to statisticians. Linear plots of response versus known concentration from each of the 9-point standard curves will be used to evaluate the linearity of the MRM measurements across the range of spiked peptide concentrations, thus providing evidence of a quantitative measurement process. Replicate analyses of the spiked plasma samples will provide estimates of assay precision (standard deviation and % CV), plus LOQ and LOD will be determined at defined signal-to-noise ratio (S/N) values. Blank runs of digested plasma with labeled peptides will provide estimates of chemical background levels in the absence of unlabeled signature peptide peaks as well as the presence of endogenous peptides in the sample (*i.e.*, CRP). Furthermore, an estimate of carryover will be determined by running a series of gradient HPLC washout runs. Finally, variation across CPTAC sites will be assessed for each of these characteristic analytical metrics.

4. Troubleshooting

This study contains a large number of samples that need to be run in a defined order and requires very reproducible peptide RTs. Problems may arise that will affect RT stability, including, but not limited to, increased column pressure, tip blockage, significant change in ambient temperature or injection of air into the system. The troubleshooting section provides some suggested routes of diagnosis and an outline of how to resolve the problem and continue with data acquisition. Importantly, operators are instructed to use the SSS and Skyline to observe any problems with retention time shifting, asymmetric peaks, low signal, and asymmetric peak shape. These problems can be caught and rectified without much down time and with minimal re-running of sample.

Procedures

1. System Suitability Sample – Column Conditioning

- a. Sample preparation
 - i. Dilute SSS from 1 pmol/ μ L to 50 fmol/ μ L
 1. Add **95 μ L** of 0.1% formic acid/3% acetonitrile to an Eppendorf tube
 2. Add 5 μ L of MichromMix (SSS) to the vial and vortex
 3. Centrifuge sample for 1 min in a benchtop centrifuge
 4. Transfer to an autosampler vial and place in autosampler
- b. Method Preparation
 - i. Prepare SSS transition list from the appropriate Skyline file
 1. 115 transitions, 22 peptides
 2. Q1, Q3 resolution = unit
 3. Dwell time = 10 msec
 4. Interscan delay (where applicable) = 5 msec (3 msec for AB Sciex 5500)
 - ii. Use the **Study 9S LC gradient** (also see Table 2b)
 1. 0-5 min, 3% B; 5-8 min, 3-15% B; 8-42 min 15-35% B; 42-45 min, 35-90% B, 45-49 min, 90% B hold; 49-50 min, 90-3% B; 50-80 min, 3% B.
Note: operators with Waters Aquity UPLC will use a slightly adjusted gradient (due to the use of Mobile phase B: 100% acetonitrile (v/v)/0.1% formic acid in water (v/v) only for Waters UPLC, see Table 3b).
 2. Flow rate: 300 nL/min
 3. Mobile phase A: 0.1% formic acid in water (v/v); Mobile phase B: 90% acetonitrile (v/v)/0.1% formic acid in water (v/v). *Note: operators with Waters Aquity UPLC will use Mobile phase B: 100% acetonitrile (v/v)/0.1% formic acid in water (v/v) only for Waters UPLC).*
 - iii. Autosampler Settings
 1. Direct Injection, full-loop injection, 1 μ L sample loop
 2. Pick up 2 μ L of sample sandwiched between mobile phase A (0.1% formic acid)
 3. See Table 1 for example
 - iv. Column heater Settings
 1. If using a column heater, please set to 35 degrees Celsius.
 - v. Recommended MS Source Conditions
 1. AB Sciex (4000 QTRAP and QTRAP 5500)
 - a. ESI voltage: 2200 \pm 200 V ; Curtain gas: 20 \pm 5; GS1: 5 \pm 5; IHT: 150 $^{\circ}$ C
 2. ThermoFisher Scientific (Vantage)
 - b. ESI voltage: 1300-1400 V; Declustering voltage: -1 V; Capillary Temperature: 210 $^{\circ}$ C
 3. Waters Xevo
 - c. Capillary voltage: 2800 V; Cone voltage: 35V; Ion source temperature: 150 $^{\circ}$ C; Cone gas flow: 15 L/Hr; Nanoflow gas flow: 0.20 bar
 4. Agilent 6410/6460

- d. ESI voltage: 1700-1800 V; gas flow: 2.5 L/min; temperature: 325°C; fragmentor voltage: 125 V
 - For all instruments:
 - Q1, Q3 resolution = unit
 - Dwell time = 10 msec
 - Interscan delay (where applicable) = 5 msec (3 msec for AB Sciex 5500)
- c. Inject sample 5 times
- d. Data Analysis
 - i. Import data files into Skyline System Suitability Sample file
 - ii. Check automatic integration of all peaks
 - 1. Manually adjust integration of peaks, if necessary
 - 2. Make sure integration start and stop is identical for all transitions of a precursor (go to “Settings”, and check “Integrate All” to enable this feature automatically)
 - 3. Document any peak that looks problematic by utilizing the “note” feature in Skyline
 - a. Excessive tailing or fronting
 - b. Drop-out of electrospray
 - c. Poor peak shape
 - d. Missing transitions
 - iii. Check that peaks pass criteria
 - 1. RT shift is isolated to first 3 injections
 - 2. Peak area CV is less than 30% for all peaks
 - 3. All peaks are detected in the last 2 sample runs
 - iv. If data do not pass criteria, troubleshoot the LC system and re-run SSS column conditioning procedure
 - v. If necessary, contact CPTAC VWG members for advice or assistance
 - 1. Susan Abbatiello, susana@broadinstitute.org
 - 2. Birgit Schilling, bschilling@buckinstitute.org

2. Retention Time Scheduling for Heavy Peptides

- a. Sample preparation
 - i. Dilute Heavy-Only IS peptide mixture from 800 fmol/μL to 80 fmol/μL
 - 1. Add 45 μL 3% acetonitrile/0.1% formic acid to an Eppendorf tube
 - 2. Add 5 μL of the 800 fmol/mL Heavy-Only IS peptide stock to the vial, vortex
 - 3. Transfer to an autosampler vial and place in autosampler
- b. Method Preparation
 - i. Prepare the Heavy-Only IS peptide transition list from the appropriate Skyline file
 - 1. 369 transitions, 123 peptides, 123 precursors
 - 2. Q1, Q3 resolution = unit
 - 3. Dwell time = 10 msec (can be lower for AB Sciex 5500 instrument, minimum is 5 msec)
 - 4. Interscan delay (where applicable) = 5 msec (3 msec for AB Sciex 5500)
 - 5. Maximum number of concurrent transitions: 60
 - 6. Export transition list as multiple methods, ignore proteins

- ii. Use the **Study 9-1 LC gradient** (also see Table 2a)
 1. 0-5 min, 3% B; 5-8 min, 3-7% B; 8-35 min 7-25% B; 35-42 min, 25-40% B; 42-45 min, 40-90% B; 45-49 min, 90% B hold; 49-50 min, 90-3% B; 50-80 min, 3% B. *Note: operators with Waters Aquity UPLC will use a slightly adjusted gradient (due to the use of Mobile phase B: 100% acetonitrile (v/v)/0.1% formic acid in water (v/v) only for Waters UPLC, see Table 3a).*
 2. Flow rate: 300 nL/min
 3. Mobile phase A: 0.1% formic acid in water (v/v); Mobile phase B: 90% acetonitrile (v/v)/0.1% formic acid in water (v/v). *Note: operators with Waters Aquity UPLC will use Mobile phase B: 100% acetonitrile (v/v)/0.1% formic acid in water (v/v) only for Waters UPLC).*
- iii. Autosampler Settings
 1. Direct Injection, full-loop injection, 1 µL sample load
 2. Pick up 2 µL of sample sandwiched between mobile phase A (0.1% formic acid)
 3. See Table 1 for example
- iv. Recommended MS Source Conditions (use as guidelines, see above, Section 1.b.iv.)
- c. Inject the sample 6 times, one time for each method created in Skyline
 - i. Name each injection with its own filename/samplename so that each injection is its own, individual file
 1. Study9-1_SiteXX_HeavyOnly_method01
 2. Study9-1_SiteXX_HeavyOnly_method02
 3. Study9-1_SiteXX_HeavyOnly_method03
 4. Study9-1_SiteXX_HeavyOnly_method04
 5. Study9-1_SiteXX_HeavyOnly_method05
 6. Study9-1_SiteXX_HeavyOnly_method06
 - ii. Follow with 1-5 injections of SSS
 1. Do not allow the system to sit idle after the 6 HeavyOnly sample runs
- d. Data Analysis
 - i. Import data files into the appropriate Skyline Platform-specific Study 9-1 file
 1. Import multiple files as one replicate
 2. Name replicate (HeavyOnly_rep1)
 - ii. Check automatic integration of all peaks
 1. Manually adjust integration of peaks, if necessary
 2. Make sure integration start and stop is identical for all transitions of a precursor (check integrate all in Skyline)
 3. Document any peak that looks problematic using the “note” feature in Skyline - Use Skyline **Custom Annotation** to note any of the above features as described below (see separate document)
 - a. Excessive tailing or fronting
 - b. Drop-out of electrospray
 - c. Poor peak shape
 - d. Missing transitions

- iii. Under Peptide Settings, set window to 4 minutes
 - 1. Export the transition list in scheduled mode – heavy only, ~325 transitions with **4 min RT window**. Double check that all 325 transitions are exported and added into the MS method appropriately.
 - 2. Run two replicates of Study 9-1 A and import into Skyline
 - 3. Ensure all peaks are detected
- iv. Under Peptide Settings, **set RT window to 2 minutes**.
- v. Under Edit->Refine, add light labeled peptides to the document
 - 1. There should now be 738 transitions, 123 peptides, 246 precursors
 - 2. View the RT plot for Scheduling
 - 3. Ensure the maximum number of concurrent transitions throughout the gradient is < 80
 - 4. If there are incidents of concurrent transitions > 80 anywhere in the gradient, contact CPTAC VWG members for assistance
- e. Export all transitions (~750) in Scheduled mode, with 2 min RT window. Double check that all 738 transitions are exported and added into the MS method appropriately.
- f. Acquire 4 replicates of sample Study9-1 A with the scheduled LC-MRM-MS method
 - i. Import data back into Skyline and observe any shift in retention times
 - ii. Use RT graph, peptide replicate view
 - iii. If peaks shift by > 30 seconds (0.5 minutes), do not proceed and contact CPTAC VWG members for assistance
 - iv. Look at data to make sure peak apexes are in or near the center of the RT window and not getting cut off during detection.
- g. Continue acquiring Study9-1 samples according to the sample chart shown below.

Important Note:

Different than in Study 7, the four replicate concentration points will be acquired in **4x singlicate curves** (with blanks and SSS in between to avoid carryover). Also note each singlicate concentration point was prepared individually at the central preparation site (Vanderbilt), so for each concentration point there will be 4 tubes supplied.

Adding Calibration Curve samples and Blinded Samples into Autosampler Vials:

Add one Calibration Curve and corresponding Blinded Samples at a time into the Autosampler. Each Calibration Curve concentration point tube contains 25 µl sample, for each singlicate curve remove the appropriate singlicate aliquot tube from the freezer, fill 12 µl into an autosampler vial, and store the remainder of that aliquot in the +4 degrees centigrade refrigerator (to avoid additional freeze-thaw cycles). Make sure to finish each singlicate curve (including initial data analysis and review of data points for quality, or for retention time drift of peaks outside the RT scheduling window) before starting the next singlicate curve. If data points need to be rerun from the current singlicate curve use remaining sample that is still in autosampler vials or if multiple reruns are needed use the saved other half of the specific singlicate aliquot that was transferred into the refrigerator as described above. For the 9-1-A sample used as blank (A-sample) at the beginning of a singlicate curve, and as blank in between blinded samples fill the entire A aliquot into an autosampler vial.

In Case of Acquisition Problems:

In case there are problems with the data acquisition or RT drift, and possible reruns that may become necessary, please keep the run number at the end of the filename to keep track of the run order. The run number should agree with the chronological order in which the samples were analyzed by LC-MRM-MS.

Blinded Samples:

Six blinded samples will be run with 1 wash run (10 fmol 6 protein mix) and 1 blank run (heavy peptides in plasma matrix, 9-1-A) in between. A total of 4 replicates of the 6 blinded samples will be acquired, one set of six in between each of the four singlicate calibration curves yielding a total of 24 blinded samples (Samples 9-1-01Blinded to 9-1-24Blinded). The wash runs use a different wash gradient **9-1-Wash** (see Table 2c below).

For Blinded Samples use Study 9-1 gradient (Table 2a) and "Study 9 scheduled MRM-MS, light and heavy" MS-method, for Wash Samples use Wash gradient (Table 2c) and "Study 9S, SSS"-adjusted MS-method, and for Blank Samples use Study 9-1 gradient (Table 2a) and "Study 9 scheduled MRM-MS, light and heavy" MS-method. Blanks (in between blinded samples) can be injected from the same autosampler vial.

Fresh Preparation of System Suitability samples (SSS) for each Singlicate Curve:

For each Singlicate Calibration Curve a fresh autosampler vial of the SSS will be prepared to prevent degradation of the sample (note in a notebook when fresh SSS are added).

Please adhere closely to sample naming as defined in the run order (with chronological runorder number at the end).

Please remember: there are 3 different gradients to be used in the below run orders:

- **Study 9S** gradient for System Suitability Samples
- **Study 9-1** gradient for the main study (all concentration curve points as well as blinded samples)
- **Study 9 wash** gradient (in between blinded samples)

Gradients are described in detail in Tables 2a-c, and 3a-c for Waters UPLC (see below)

First singlicate curve plus blinded samples (01-06blinded)

| Run Number | Number of injections | Sample Description | Filename | Method | Notes |
|------------|----------------------|------------------------------------|---|--|--|
| 1 | 1 | SSS, Michrom Mix | 9-1 SiteX SSS run 001 | Study 9S (SSS) | Column Conditioning |
| 2 | 1 | SSS, Michrom Mix | 9-1 SiteX SSS run 002 | Study 9S (SSS) | |
| 3 | 1 | SSS, Michrom Mix | 9-1 SiteX SSS run 003 | Study 9S (SSS) | |
| 4 | 1 | SSS, Michrom Mix | 9-1 SiteX SSS run 004 | Study 9S (SSS) | |
| 5 | 1 | SSS, Michrom Mix | 9-1 SiteX SSS run 005 | Study 9S (SSS) | |
| 6 | 1 | Study9-1 Heavy Only | 9-1 SiteX heavy method01 run 006 | Study 9 method 01 | For scheduling |
| 7 | 1 | Study9-1 Heavy Only | 9-1 SiteX heavy method02 run 007 | Study 9 method 02 | |
| 8 | 1 | Study9-1 Heavy Only | 9-1 SiteX heavy method03 run 008 | Study 9 method 03 | |
| 9 | 1 | Study9-1 Heavy Only | 9-1 SiteX heavy method04 run 009 | Study 9 method 04 | |
| 10 | 1 | Study9-1 Heavy Only | 9-1 SiteX heavy method05 run 010 | Study 9 method 05 | |
| 11 | 1 | Study9-1 Heavy Only | 9-1 SiteX heavy method06 run 011 | Study 9 method 06 | |
| 12-16 | 5 | SSS, Michrom Mix | 9-1_SiteX_SSS_run_012 9-1_SiteX_SSS_run_013 9-1_SiteX_SSS_run_014 9-1_SiteX_SSS_run_015 9-1_SiteX_SSS_run_016 | Study 9S (SSS) | At least 3 SSS samples and if necessary place holder |
| 17 | 1 | Study9-1 A, IS peptides in plasma | 9-1_SiteX_A_sMRM_run_017 | Study 9 scheduled MRM-MS, heavy only | Scheduled run, 4 min RT window |
| 18 | 1 | Study9-1 A, IS peptides in plasma | 9-1_SiteX_A_sMRM_run_018 | Study 9 scheduled MRM-MS, heavy only | |
| 19 | 1 | Study9-1 A1, IS peptides in plasma | 9-1_SiteX_A1_CalCurve_run_019 | Study 9 scheduled MRM-MS, light and heavy | Scheduled run, 2 min RT window, Calibration Curve |
| 20 | 1 | Study9-1 B1 | 9-1_SiteX_B1_CalCurve_run_020 | Study 9 scheduled MRM-MS, light and heavy | |
| 21 | 1 | Study9-1 C1 | 9-1_SiteX_C1_CalCurve_run_021 | Study 9 scheduled MRM-MS, light and heavy, "Study9 sMRM" | CalCurve |
| 22 | 1 | Study9-1 D1 | 9-1 SiteX D1 CalCurve run 022 | Study9 sMRM | CalCurve |
| 23 | 1 | SSS | 9-1_SiteX_SSS_run_023 | Study 9S (SSS) | System Suitability |
| 24 | 1 | Study9-1 E1 | 9-1_SiteX_E1_CalCurve_run_024 | Study9 sMRM | CalCurve |

| | | | | | |
|----|---|-----------------------------------|--------------------------------|---|--------------------|
| 25 | 1 | Study9-1 F1 | 9-1_SiteX_F1_CalCurve_run_025 | Study9 sMRM | CalCurve |
| 26 | 1 | Study9-1 G1 | 9-1_SiteX_G1_CalCurve_run_026 | Study9 sMRM | CalCurve |
| 27 | 1 | Study9-1 H1 | 9-1_SiteX_H1_CalCurve_run_027 | Study9 sMRM | CalCurve |
| 28 | 1 | Study9-1 QC1 | 9-1_SiteX_QC1_CalCurve_run_028 | Study9 sMRM | QC |
| 29 | 1 | Study9-1 I1 | 9-1_SiteX_I1_CalCurve_run_029 | Study9 sMRM | CalCurve |
| 30 | 1 | Study9-1 J1 | 9-1_SiteX_J1_CalCurve_run_030 | Study9 sMRM | CalCurve |
| 31 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_31 | Study 9 wash | |
| 32 | 1 | SSS | 9-1_SiteX_SSS_run_032 | Study 9S (SSS) | System suitability |
| 33 | 1 | SSS | 9-1_SiteX_SSS_run_033 | Study 9S (SSS) | System suitability |
| 34 | 1 | SSS | 9-1_SiteX_SSS_run_034 | Study 9S (SSS) | System suitability |
| 35 | 1 | Study9-1 A, IS peptides in plasma | 9-1_SiteX_A_blank_run_035 | Study 9 scheduled MRM-MS, light and heavy | |
| 36 | 1 | Study9-1 01blinded-Sample | 9-1_SiteX_01blinded_run_036 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 37 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_037 | Study 9 wash | |
| 38 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_038 | Study 9 scheduled MRM-MS, light and heavy | |
| 39 | 1 | Study9-1 02blinded | 9-1_SiteX_02blinded_run_039 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 40 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_040 | Study 9 wash | |
| 41 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_041 | Study 9 scheduled MRM-MS, light and heavy | |
| 42 | 1 | SSS | 9-1_SiteX_SSS_run_042 | Study 9S (SSS) | System suitability |
| 43 | 1 | Study9-1 03blinded | 9-1_SiteX_03blinded_run_043 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 44 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_044 | Study 9 wash | |

| | | | | | |
|----|---|------------------------------|--|---|--------------------|
| 45 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_045 | Study 9 scheduled MRM-MS, light and heavy | |
| 46 | 1 | Study9-1 04blinded | 9-1_SiteX_04blinded_run_046 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 47 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_047 | Study 9 wash | |
| 48 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_048 | Study 9 scheduled MRM-MS, light and heavy | |
| 49 | 1 | SSS | 9-1_SiteX_SSS_run_049 | Study 9S (SSS) | System suitability |
| 50 | 1 | Study9-1 05blinded | 9-1_SiteX_05blinded_run_050 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 51 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_051 | Study 9 wash | |
| 52 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_052 | Study 9 scheduled MRM-MS, light and heavy | |
| 53 | 1 | Study9-1 06blinded | 9-1_SiteX_06blinded_run_053 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 54 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_054 | Study 9 wash | |
| 55 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_055 | Study 9 scheduled MRM-MS, light and heavy | |
| | | | Run second, third and fourth singlicate curves (see below) | | |

Second singlicate curve plus blinded samples (07-12blinded)

| | | | | | |
|-------|---|------------------------------------|---|--|---|
| 56-58 | 3 | SSS, Michrom Mix | 9-1_SiteX_SSS_run_056 9-1_SiteX_SSS_run_057 9-1_SiteX_SSS_run_058 | Study 9S (SSS) | 3 SSS samples |
| 59 | 1 | Study9-1 A2, IS peptides in plasma | 9-1_SiteX_A2_CalCurve_run_059 | Study 9 scheduled MRM-MS, light and heavy | Scheduled run, 2 min RT window, Calibration Curve |
| 60 | 1 | Study9-1 B2 | 9-1_SiteX_B2_CalCurve_run_060 | Study 9 scheduled MRM-MS, light and heavy | |
| 61 | 1 | Study9-1 C2 | 9-1_SiteX_C2_CalCurve_run_061 | Study 9 scheduled MRM-MS, light and heavy, "Study9 sMRM" | CalCurve |
| 62 | 1 | Study9-1 D2 | 9-1_SiteX_D2_CalCurve_run_062 | Study9 sMRM | CalCurve |
| 63 | 1 | SSS | 9-1_SiteX_SSS_run_063 | Study 9S (SSS) | System Suitability |
| 64 | 1 | Study9-1 E2 | 9-1_SiteX_E2_CalCurve_run_064 | Study9 sMRM | CalCurve |
| 65 | 1 | Study9-1 F2 | 9-1_SiteX_F2_CalCurve_run_065 | Study9 sMRM | CalCurve |
| 66 | 1 | Study9-1 G2 | 9-1_SiteX_G2_CalCurve_run_066 | Study9 sMRM | CalCurve |
| 67 | 1 | Study9-1 H2 | 9-1_SiteX_H2_CalCurve_run_067 | Study9 sMRM | CalCurve |
| 68 | 1 | Study9-1 QC2 | 9-1_SiteX_QC2_CalCurve_run_068 | Study9 sMRM | QC |
| 69 | 1 | Study9-1 I2 | 9-1_SiteX_I2_CalCurve_run_069 | Study9 sMRM | CalCurve |
| 70 | 1 | Study9-1 J2 | 9-1_SiteX_J2_CalCurve_run_070 | Study9 sMRM | CalCurve |
| 71 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_071 | Study 9 wash | |
| 72 | 1 | SSS | 9-1_SiteX_SSS_run_072 | Study 9S (SSS) | System suitability |
| 73 | 1 | SSS | 9-1_SiteX_SSS_run_073 | Study 9S (SSS) | System suitability |
| 74 | 1 | SSS | 9-1_SiteX_SSS_run_074 | Study 9S (SSS) | System suitability |
| 75 | 1 | Study9-1 A, IS peptides in plasma | 9-1_SiteX_A_blank_run_075 | Study 9 scheduled MRM-MS, light and heavy | |
| 76 | 1 | Study9-1 07blinded-Sample | 9-1_SiteX_07blinded_run_076 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 77 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_077 | Study 9 wash | |

| | | | | | |
|----|---|------------------------------|-----------------------------|---|--------------------|
| 78 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_078 | Study 9 scheduled MRM-MS, light and heavy | |
| 79 | 1 | Study9-1 08blinded | 9-1_SiteX_08blinded_run_079 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 80 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_080 | Study 9 wash | |
| 81 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_081 | Study 9 scheduled MRM-MS, light and heavy | |
| 82 | 1 | SSS | 9-1_SiteX_SSS_run_082 | Study 9S (SSS) | System suitability |
| 83 | 1 | Study9-1 09blinded | 9-1_SiteX_09blinded_run_083 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 84 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_084 | Study 9 wash | |
| 85 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_085 | Study 9 scheduled MRM-MS, light and heavy | |
| 86 | 1 | Study9-1 10blinded | 9-1_SiteX_10blinded_run_086 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 87 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_087 | Study 9 wash | |
| 88 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_088 | Study 9 scheduled MRM-MS, light and heavy | |
| 89 | 1 | SSS | 9-1_SiteX_SSS_run_089 | Study 9S (SSS) | System suitability |
| 90 | 1 | Study9-1 11blinded | 9-1_SiteX_11blinded_run_090 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 91 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_091 | Study 9 wash | |
| 92 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_092 | Study 9 scheduled MRM-MS, light and heavy | |

| | | | | | |
|----|---|------------------------------|--|---|---------|
| 93 | 1 | Study9-1 12blinded | 9-1_SiteX_12blinded_run_093 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 94 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_094 | Study 9 wash | |
| 95 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_095 | Study 9 scheduled MRM-MS, light and heavy | |
| | | | Run third and fourth singlicate curves (see below) | | |

Third singlicate curve plus blinded samples (13-18blinded)

| | | | | | |
|-------|---|------------------------------------|---|--|---|
| 96-98 | 3 | SSS, Michrom Mix | 9-1_SiteX_SSS_run_096 9-1_SiteX_SSS_run_097 9-1_SiteX_SSS_run_098 | Study 9S (SSS) | 3 SSS samples |
| 99 | 1 | Study9-1 A3, IS peptides in plasma | 9-1_SiteX_A3_CalCurve_run_099 | Study 9 scheduled MRM-MS, light and heavy | Scheduled run, 2 min RT window, Calibration Curve |
| 100 | 1 | Study9-1 B3 | 9-1_SiteX_B3_CalCurve_run_100 | Study 9 scheduled MRM-MS, light and heavy | |
| 101 | 1 | Study9-1 C3 | 9-1_SiteX_C3_CalCurve_run_101 | Study 9 scheduled MRM-MS, light and heavy, "Study9 sMRM" | CalCurve |
| 102 | 1 | Study9-1 D3 | 9-1_SiteX_D3_CalCurve_run_102 | Study9 sMRM | CalCurve |
| 103 | 1 | SSS | 9-1_SiteX_SSS_run_103 | Study 9S (SSS) | System Suitability |
| 104 | 1 | Study9-1 E3 | 9-1_SiteX_E3_CalCurve_run_104 | Study9 sMRM | CalCurve |
| 105 | 1 | Study9-1 F3 | 9-1_SiteX_F3_CalCurve_run_105 | Study9 sMRM | CalCurve |
| 106 | 1 | Study9-1 G3 | 9-1_SiteX_G3_CalCurve_run_106 | Study9 sMRM | CalCurve |
| 107 | 1 | Study9-1 H3 | 9-1_SiteX_H3_CalCurve_run_107 | Study9 sMRM | CalCurve |
| 108 | 1 | Study9-1 QC3 | 9-1_SiteX_QC3_CalCurve_run_108 | Study9 sMRM | QC |
| 109 | 1 | Study9-1 I3 | 9-1_SiteX_I3_CalCurve_run_109 | Study9 sMRM | CalCurve |
| 110 | 1 | Study9-1 J3 | 9-1_SiteX_J3_CalCurve_run_110 | Study9 sMRM | CalCurve |

| | | | | | |
|-----|---|-----------------------------------|-----------------------------|---|--------------------|
| 111 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_111 | Study 9 wash | |
| 112 | 1 | SSS | 9-1_SiteX_SSS_run_112 | Study 9S (SSS) | System suitability |
| 113 | 1 | SSS | 9-1_SiteX_SSS_run_113 | Study 9S (SSS) | System suitability |
| 114 | 1 | SSS | 9-1_SiteX_SSS_run_114 | Study 9S (SSS) | System suitability |
| 115 | 1 | Study9-1 A, IS peptides in plasma | 9-1_SiteX_A_blank_run_115 | Study 9 scheduled MRM-MS, light and heavy | |
| 116 | 1 | Study9-1 13blinded-Sample | 9-1_SiteX_13blinded_run_116 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 117 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_117 | Study 9 wash | |
| 118 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_118 | Study 9 scheduled MRM-MS, light and heavy | |
| 119 | 1 | Study9-1 14blinded | 9-1_SiteX_14blinded_run_119 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 120 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_120 | Study 9 wash | |
| 121 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_121 | Study 9 scheduled MRM-MS, light and heavy | |
| 122 | 1 | SSS | 9-1_SiteX_SSS_run_122 | Study 9S (SSS) | System suitability |
| 123 | 1 | Study9-1 15blinded | 9-1_SiteX_15blinded_run_123 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 124 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_124 | Study 9 wash | |
| 125 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_125 | Study 9 scheduled MRM-MS, light and heavy | |
| 126 | 1 | Study9-1 16blinded | 9-1_SiteX_16blinded_run_126 | Study 9 scheduled MRM-MS, light | Blinded |

| | | | | | |
|-----|---|------------------------------|---|---|--------------------|
| | | | | and heavy | |
| 127 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_127 | Study 9 wash | |
| 128 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_128 | Study 9 scheduled MRM-MS, light and heavy | |
| 129 | 1 | SSS | 9-1_SiteX_SSS_run_129 | Study 9S (SSS) | System suitability |
| 130 | 1 | Study9-1 17blinded | 9-1_SiteX_17blinded_run_130 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 131 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_131 | Study 9 wash | |
| 132 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_132 | Study 9 scheduled MRM-MS, light and heavy | |
| 133 | 1 | Study9-1 18blinded | 9-1_SiteX_18blinded_run_133 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 134 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_134 | Study 9 wash | |
| 135 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_135 | Study 9 scheduled MRM-MS, light and heavy | |
| | | | Run fourth singlicate curve (see below) | | |

Fourth singlicate curve plus blinded samples (19-24blinded)

| | | | | | |
|---------|---|------------------------------------|---|---|---|
| 136-138 | 3 | SSS, Michrom Mix | 9-1_SiteX_SSS_run_136 9-1_SiteX_SSS_run_137 9-1_SiteX_SSS_run_138 | Study 9S (SSS) | 3 SSS samples |
| 139 | 1 | Study9-1 A4, IS peptides in plasma | 9-1_SiteX_A4_CalCurve_run_139 | Study 9 scheduled MRM-MS, light and heavy | Scheduled run, 2 min RT window, Calibration Curve |
| 140 | 1 | Study9-1 B4 | 9-1_SiteX_B4_CalCurve_run_140 | Study 9 scheduled MRM-MS, light and heavy | |

| | | | | | |
|-----|---|-----------------------------------|--------------------------------|--|--------------------|
| 141 | 1 | Study9-1 C4 | 9-1_SiteX_C4_CalCurve_run_141 | Study 9 scheduled MRM-MS, light and heavy, “Study9 sMRM” | CalCurve |
| 142 | 1 | Study9-1 D4 | 9-1_SiteX_D4_CalCurve_run_142 | Study9 sMRM | CalCurve |
| 143 | 1 | SSS | 9-1_SiteX_SSS_run_143 | Study 9S (SSS) | System Suitability |
| 144 | 1 | Study9-1 E4 | 9-1_SiteX_E4_CalCurve_run_144 | Study9 sMRM | CalCurve |
| 145 | 1 | Study9-1 F4 | 9-1_SiteX_F4_CalCurve_run_145 | Study9 sMRM | CalCurve |
| 146 | 1 | Study9-1 G4 | 9-1_SiteX_G4_CalCurve_run_146 | Study9 sMRM | CalCurve |
| 147 | 1 | Study9-1 H4 | 9-1_SiteX_H4_CalCurve_run_147 | Study9 sMRM | CalCurve |
| 148 | 1 | Study9-1 QC4 | 9-1_SiteX_QC4_CalCurve_run_148 | Study9 sMRM | QC |
| 149 | 1 | Study9-1 I4 | 9-1_SiteX_I4_CalCurve_run_149 | Study9 sMRM | CalCurve |
| 150 | 1 | Study9-1 J4 | 9-1_SiteX_J4_CalCurve_run_150 | Study9 sMRM | CalCurve |
| 151 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_151 | Study 9 wash | |
| 152 | 1 | SSS | 9-1_SiteX_SSS_run_152 | Study 9S (SSS) | System suitability |
| 153 | 1 | SSS | 9-1_SiteX_SSS_run_153 | Study 9S (SSS) | System suitability |
| 154 | 1 | SSS | 9-1_SiteX_SSS_run_154 | Study 9S (SSS) | System suitability |
| 155 | 1 | Study9-1 A, IS peptides in plasma | 9-1_SiteX_A_blank_run_155 | Study 9 scheduled MRM-MS, light and heavy | |
| 156 | 1 | Study9-1 19blinded-Sample | 9-1_SiteX_19blinded_run_156 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 157 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_157 | Study 9 wash | |
| 158 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_158 | Study 9 scheduled MRM-MS, light and heavy | |
| 159 | 1 | Study9-1 20blinded | 9-1_SiteX_20blinded_run_159 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 160 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_160 | Study 9 wash | |
| 161 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_161 | Study 9 scheduled | |

| | | | | | |
|-----|---|------------------------------|-----------------------------|---|--------------------|
| | | | | MRM-MS, light and heavy | |
| 162 | 1 | SSS | 9-1_SiteX_SSS_run_162 | Study 9S (SSS) | System suitability |
| 163 | 1 | Study9-1 21blinded | 9-1_SiteX_21blinded_run_163 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 164 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_164 | Study 9 wash | |
| 165 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_165 | Study 9 scheduled MRM-MS, light and heavy | |
| 166 | 1 | Study9-1 22blinded | 9-1_SiteX_22blinded_run_166 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 167 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_167 | Study 9 wash | |
| 168 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_168 | Study 9 scheduled MRM-MS, light and heavy | |
| 169 | 1 | SSS | 9-1_SiteX_SSS_run_169 | Study 9S (SSS) | System suitability |
| 170 | 1 | Study9-1 23blinded | 9-1_SiteX_23blinded_run_170 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 171 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_171 | Study 9 wash | |
| 172 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_172 | Study 9 scheduled MRM-MS, light and heavy | |
| 173 | 1 | Study9-1 24blinded | 9-1_SiteX_24blinded_run_173 | Study 9 scheduled MRM-MS, light and heavy | Blinded |
| 174 | 1 | Study9-1 wash | 9-1_SiteX_wash_run_174 | Study 9 wash | |
| 175 | 1 | IS peptides in plasma, blank | 9-1_SiteX_A_blank_run_175 | Study 9 scheduled MRM-MS, light and heavy | |
| | | | Study 9-1finished !! | | |

1. Data Analysis

- a. All data should be analyzed in “real time,” or as close to its acquisition time as possible
- b. Data analysis must be conducted through Skyline using templates provided for each instrument platform
- c. Once data are imported, check integration of all peaks
 - i. Integrate transitions with the same start and stop time: Under “Settings” in Skyline, enable “Integrate All”
 - ii. Integrate light and heavy peptides with the same start and stop time
 - iii. If an interference is present, document the transition-level note for that transition
- d. Export data for submission to VWG statisticians using the report template provided
- e. Under Settings, click on “custom annotation” and make sure that all settings are selected.
 - i. If a transition, precursor, or peptide is deemed unusable because of a bad injection, loss of electrospray or another explainable reason and the sample is re-run, please click the “do not use” box for the sample that is to be excluded from subsequent analysis for LOD/LOQ.

Even if a specific replicate needs to be rerun, the original and new (rerun) replicate should both be imported into Skyline, and "custom annotation" capabilities of Skyline should be used to annotate what data points should be used for statistical analysis.

4. Trouble Shooting Section

- a. **Symptom: No data (missing peaks) or poor quality peaks for Study 9 peptides**
 - i. Look at pressure trace to make sure sample (not air) was injected onto column.
 - ii. Open up the last SSS acquired before file with missing peaks
 1. In Skyline or in vendor specific software
 - iii. Are 22 peaks present?
 1. Yes: go on to point iv.
 2. No: identify which peaks are missing
 - a. Check pressure trace to make sure air was not injected onto column and that column pressure isn't too high
 - b. Is there enough sample in all vials (sample vials and reagent vials)?
 - iv. Are peak shapes symmetric and well defined as opposed to mis-shapen and jagged?
 1. Yes: go on to point v.
 2. No: check ESI tip of column, clean if necessary.
 - a. Also check pressure trace to make sure pressure is not too high.
 - v. Is column pressure too high?
 1. Compare column pressure of bad runs to previous runs in the sample list
 - a. If pressure is too high, remove column and clean tip
 - i. Wipe tip with gloved hand, wet with water and wipe again. Check tip under microscope if possible.
 - ii. If this does not work, check for blockages in transfer lines closest to column and work backwards to LC.
 - iii. If this does not work, consider replacing column.

- vi. Is column pressure too low?
 - 1. If pressure is too low, check all connections between column and LC.
- vii. Are flow rates calibrated correctly? (potentially check or recalibrate flow rates)
- viii. Are samples being picked up properly from the autosampler vials?

b. Symptom: Retention Time Drift

- i. Observe RT shift of the SSS samples in a single Skyline document
 - 1. Does RT shift in one direction or randomly?
 - ii. If in one direction, are the peaks migrating out of the 2 min RT window for the Study 9 method? Solution: Re-scheduling because of shifting retention times
 - 1. Re-run Study 9-1_A with the “heavy only” scheduled method that has a 4 minute RT window.
 - 2. Import data into an empty Skyline template for Study 9 peptides and export a new scheduled method based on the RT data for light and heavy peptides (738 transitions).
 - 3. Keep close eye on RT shifts.
 - 4. Consider running SSS runs at points when you would not be able to check the RT shift and if peaks are in danger of shifting out of RT window (like overnight).
- c. If any concentration points must be re-run, it is important to follow this procedure:**
- i. Run SSS sample to ensure system is working properly. Look for smooth, symmetrical peaks with stable retention times.
 - ii. Inject Sample Study9-1_A before re-running a calibration point to make sure there is little carryover.
 - iii. Follow with another SSS run so that you have time to analyze the data before running the next concentration point.

HPLC Chromatography Conditions for Studies 9-1:

Individual CPTAC sites are expected to implement these HPLC conditions for the duration of the study. [Packed New Objective HPLC Columns are provided.](#)

- HPLC-plumbing: DIRECT INJECTION
- Sample Loop: 1 µL sample loop
- Columns for AB and Thermo HPLC systems: (New Objective custom-packed as previously shipped to all sites): PicoFrit 75 µm ID / 10 µm tip (ReproSil-Pur C18-AQ, 3µm, 120Å, length of 12 cm)
- Columns for Waters UPLC systems: a Waters specific column will be used (1.7µm BEH130 C18, 75um X 150mm column)
- Columns for Agilent ChipCube systems: NanoChip
- Mobile phases: (A) 0.1% Formic acid (v/v); (B) 90% Acetonitrile / 0.1% Formic acid (v/v)
- Flow rate: 300 nL/min

- Injection volume: 1 µL on column (full loop injection with overloading the loop with 2 µL)
- Loop for Eksigent LC: 1 µL PEEKsil loop provided in sample kit, 100 µm ID PEEKsil.
- Injection Amount: ~0.5 µg total protein on-column
- Gradient: for details see below, **NOTE: Study 9S system suitability has a slightly different gradient than the new Study 9-1 gradient. Please use the appropriate gradient for these different applications.**

Table 1: Eksigent/Tempo Notes: Autosampler set-up, direct injection mode (representative example)

| <u>Autosampler Program with Standard Injection*</u> | | |
|--|-------------|-----------------------------------|
| # | Function | Command |
| 1 | Output | 1-Off |
| 2 | Output | 2-Off |
| 3 | Valve | Injector Load |
| 4** | Aspirate | 10 µL Reagent-1 Speed:1 Height:5 |
| 5 | Aspirate | 2 µL Sample Speed:1 Height:2 |
| 6** | Aspirate | 2.3 µL Reagent-1 Speed:1 Height:5 |
| 7 | Output | 2-On |
| 8 | Valve | Injector Inject |
| 9 | Dispense | 14.3 µL Waste Speed:5 Height:0 |
| 10 ** | Needle Wash | 200 µL |
| 11 | End | |

NOTE: Different gradients are used for Study 9-1 calibration curve (Study 9-1 gradient) vs. system suitability samples (Study 9S gradient) vs. Wash gradient used in between blinded samples (Study 9 Wash gradient).

Table 2a. 9-1 HPLC gradient for all analyses of 123 target synthetic peptides.

| Time (min) | %A | %B | Flow Rate (nl/min) |
|------------|----|----|--------------------|
| 0 | 97 | 3 | 300 |
| 5 | 97 | 3 | 300 |
| 8 | 93 | 7 | 300 |
| 35 | 75 | 25 | 300 |
| 42 | 60 | 40 | 300 |
| 45 | 10 | 90 | 300 |
| 49 | 10 | 90 | 300 |
| 50 | 97 | 3 | 300 |
| 80 | 97 | 3 | 300 |

Mobile phase A: 0.1% formic acid (v/v) in water

Mobile phase B: 90% acetonitrile/0.1% formic acid (v/v) in water

Table 2b. 9S HPLC gradient to be used for System Suitability Sample (Michrom 6 protein mix).

| Time (min) | %A | %B | Flow Rate (nl/min) |
|------------|----|----|--------------------|
| 0 | 97 | 3 | 300 |
| 5 | 97 | 3 | 300 |
| 8 | 85 | 15 | 300 |
| 42 | 65 | 35 | 300 |
| 45 | 10 | 90 | 300 |
| 49 | 10 | 90 | 300 |
| 50 | 97 | 3 | 300 |
| 80 | 97 | 3 | 300 |

Mobile phase A: 0.1% formic acid (v/v) in water

Mobile phase B: 90% acetonitrile/0.1% formic acid (v/v) in water

Table 2c. 9-Wash HPLC gradient to be used for washing the column in between Blinded Samples.

| Time (min) | %A | %B | Flow Rate (nl/min) |
|------------|----|----|--------------------|
| 0 | 97 | 3 | 300 |
| 15 | 10 | 90 | 300 |
| 16 | 97 | 3 | 300 |
| 31 | 10 | 90 | 300 |
| 32 | 97 | 3 | 300 |
| 62 | 97 | 3 | 300 |

Mobile phase A: 0.1% formic acid (v/v) in water

Mobile phase B: 90% acetonitrile/0.1% formic acid (v/v) in water

NOTE: Total peptide amount greater than 1 µg injected onto nanoLC columns can result in poor chromatographic peak shape and poor reproducibility from run to run. The MARS-14 depleted plasma samples (Samples 9-1-A through 9-1-J, and 9-1-Blank) have been diluted such that a 1 µL injection results in approximately 0.5 µg of total protein on-column. Therefore, they should be analyzed without any additional dilution.

For Waters operators only:

Operators with Waters Aquity UPLC will use a slightly adjusted gradient, due to the use of Mobile phase B: 100% acetonitrile (v/v)/0.1% formic acid in water (v/v) only for Waters UPLC.

Adjusted Table 3a., for Waters UPLC only. Waters 9-1 HPLC gradient for all analyses of 123 target synthetic peptides.

| Time (min) | %A | %B | Flow Rate (nl/min) |
|------------|------|------|--------------------|
| 0 | 97.3 | 2.7 | 300 |
| 5 | 97.3 | 2.7 | 300 |
| 8 | 93.7 | 6.3 | 300 |
| 35 | 77.5 | 22.5 | 300 |
| 42 | 64 | 36 | 300 |
| 45 | 19 | 81 | 300 |
| 49 | 19 | 81 | 300 |
| 50 | 97.3 | 2.7 | 300 |
| 80 | 97.3 | 2.7 | 300 |

Mobile phase A: 0.1% formic acid (v/v) in water

Mobile phase B: 100% acetonitrile/0.1% formic acid (v/v) in water

Adjusted Table 3b., for Waters UPLC only. Waters 9S HPLC gradient to be used for System Suitability Sample (Michrom 6 protein mix).

| Time (min) | %A | %B | Flow Rate (nl/min) |
|------------|------|------|--------------------|
| 0 | 97.3 | 2.7 | 300 |
| 5 | 97.3 | 2.7 | 300 |
| 8 | 86.5 | 13.5 | 300 |
| 42 | 68.5 | 31.5 | 300 |
| 45 | 19 | 81 | 300 |
| 49 | 19 | 81 | 300 |
| 50 | 97.3 | 2.7 | 300 |
| 80 | 97.3 | 2.7 | 300 |

Mobile phase A: 0.1% formic acid (v/v) in water

Mobile phase B: 100% acetonitrile/0.1% formic acid (v/v) in water

Adjusted Table 3c. for Waters UPLC only. Waters 9-Wash HPLC gradient to be used for washing the column in between Blinded Samples.

| Time (min) | %A | %B | Flow Rate (nl/min) |
|------------|------|-----|--------------------|
| 0 | 97.3 | 2.7 | 300 |
| 15 | 19 | 81 | 300 |
| 16 | 97.3 | 2.7 | 300 |
| 31 | 19 | 81 | 300 |
| 32 | 97.3 | 2.7 | 300 |
| 62 | 97.3 | 2.7 | 300 |

Mobile phase A: 0.1% formic acid (v/v) in water

Mobile phase B: 100% acetonitrile/0.1% formic acid (v/v) in water